

## Terrific Broth or TB

Cat No: 40041, 40042, 40043, and 40044

Storage: room temperature for six months

### Compositions:

1.2% Peptone, 2.4% Yeast extract, 72 mM  $K_2HPO_4$ , 17 mM  $KH_2PO_4$ , and 0.4% Glycerol

### Aseptic operation:

It is best to use all of the medium once it is opened.

Wash both hands and clean the working surface by 70% ethanol.

Do not cough, sneeze, or breath into the medium.

Hold breath when opening the medium and quickly close it after transferring the medium out.

Avoid opening the medium in places such as bacterial culture rooms.

Use a biological hood if contamination persists.

### Warranty:

We guarantee the sterilization and composition of our Terrific Broth. If the product fails to meet its specs, we will promptly replace the product. All other warranties are expressly disclaimed, including but not limited to, the implied warranties of merchantability and fitness for a particular purpose. We shall not be liable for any consequential damages.